

Validation of the slaughterhouse porcine heart model for ex-situ heart perfusion studies

Selma E Kaffka genaamd Dengler,¹  Mudit Mishra,^{1,2}
Sjoerd van Tuijl,³ Saskia CA de Jager,² Joost PG Sluijter,²
Pieter A Doevendans^{4,5} and Niels P van der Kaaij¹

Perfusion
2023, Vol. 0(0) 1–9
© The Author(s) 2023



Article reuse guidelines:
sagepub.com/journals-permissions
DOI: 10.1177/02676591231152718
journals.sagepub.com/home/prf



Abstract

Introduction: To validate slaughterhouse hearts for *ex-situ* heart perfusion studies, we compared cold oxygenated machine perfusion in less expensive porcine slaughterhouse hearts ($N = 7$) to porcine hearts that are harvested following the golden standard in laboratory animals ($N = 6$).

Methods: All hearts received modified St Thomas 2 crystalloid cardioplegia prior to 4 hours of cold oxygenated machine perfusion. Hearts were perfused with homemade modified Steen heart solution with a perfusion pressure of 20–25 mmHg to achieve a coronary flow between 100–200 mL/min. Reperfusion and testing was performed for 4 hours on a normothermic, oxygenated diluted whole blood loaded heart model. Survival was defined by a cardiac output above 3 L with a mean aortic pressure above 60 mmHg.

Results: Both groups showed 100% functional survival, with laboratory hearts displaying superior cardiac function. Both groups showed similar decline in function over time.

Conclusion: We conclude that the slaughterhouse heart can be used as an alternative to laboratory hearts and provides a cost-effective method for future *ex-situ* heart perfusion studies.

Keywords

cold oxygenated machine perfusion, slaughterhouse porcine hearts, *ex-situ* heart perfusion, working heart model

Introduction

With the reintroduction of donation after circulatory death (DCD) heart donation,^{1,2} both clinical and experimental research in *ex-situ* heart perfusion has intensified. A downside of experimental heart perfusion studies is the use of hearts that are harvested from laboratory (LAB) animals, of which the porcine heart is most commonly used. With the current focus on refinement, reduction and replacement in animal studies (the 3R's),³ alternatives to sacrifice laboratory animals are desired. The slaughterhouse (SLH) porcine heart can be used as an alternative source for *ex-situ* heart perfusion studies as their access is unlimited without sacrificing animals solely for research goals and thereby reducing experimental costs. The SLH-hearts have been used previously for similar studies.⁴ Although an alternative, the use of slaughterhouse porcine hearts comes at a cost which is seen in terms of higher variability due to different farms of origin, stress due to transport and loading and a fattening diet.^{5,6} Furthermore,

animals are 6 months old and around 120 kg as they are bred for human consumption, which also results in more hypertrophic hearts compared to experimental animals, and might lead to a larger variation in results obtained. Moreover, the harvesting of the SLH heart itself carries important implications. As regulations for the slaughter process require the sedation of animal

¹Department of Cardiothoracic Surgery, University Medical Center Utrecht, Utrecht, The Netherlands

²Laboratory of Experimental Cardiology, Department of Cardiology, University Medical Center Utrecht, University Utrecht, Utrecht, The Netherlands

³LifeTec Group B.V., Kennedyplein 10-11, The Netherlands

⁴Department of Cardiology, University Medical Center Utrecht, University Utrecht, Utrecht, The Netherlands

⁵Netherlands Heart Institute, Moreelsepark 1, The Netherlands

Corresponding author:

Selma E Kaffka genaamd Dengler, University Medical Center Utrecht, Postbox 85500, 3508 GA Utrecht, Heidelberglaan 100, The Netherlands.
Email: s.e.kaffkagenaamddengler@umcutrecht.nl

prior to exsanguination,⁷ the pig is electrically stunned and harvesting of the heart can only start after exsanguination of the animal. This results in an extended period of warm ischemia with potentially more severe ischemia reperfusion injury (IRI) upon reperfusion, which is more significant in hypertrophic hearts.⁸ These factors together may impact myocardial function, but the exact effect in the slaughterhouse setting is unclear.

The harvesting procedure of the heart in LAB animals is in great contrast to SLH animals, since in LAB animals the heart can be harvested in the most optimal way. Animals undergo anesthesia prior to sternotomy followed by a controlled cardioplegic arrest of the heart. This procedure minimizes warm ischemia and prevents right ventricular distension and is known as the golden standard of heart harvesting.⁹

In the current study, we compare cold oxygenated machine-perfused porcine slaughterhouse hearts to cold oxygenated machine-perfused porcine hearts that are harvested following the golden standard in laboratory animals in order to validate the slaughterhouse model for heart perfusion studies. The aim was to show that slaughterhouse hearts approach functional outcome and survival.

Material and methods

Slaughterhouse animals (SLH)

7 hearts were harvested from Dutch Landrace Hybrid pigs slaughtered for human consumption. The protocols were consistent with the EC regulations 1069/2009 regarding slaughterhouse animal material for diagnosis and research as supervised by the Dutch Government (Dutch Ministry of Agriculture, Nature and Food Quality) and were approved by the associated legal authorities of animal welfare (Food and Consumer Product Safety Authority).

Heart harvesting. Animals were electrically stunned, hung, exsanguinated, and declared dead prior to harvesting the heart. Death was determined if no reflex was seen with eye stimulus after draining at least 2 L of blood. No heparin was administered to the animal due to applicable regulations. Blood (3–3.5 L) was collected in a bag with 5000 IU/L heparin as anticoagulant. After parasternal incision, the heart and lungs were harvested en-bloc by transecting the trachea, the aorta distally from the arch vessels and the caval veins. Topical cooling was applied while the aorta was cannulated and modified St Thomas 2 crystalloid cardioplegia was administered at a pressure of 80–100 mmHg.⁴ Warm ischemic time recorded ranged between the stunning of the animal to the delivery of cardioplegia.

Laboratory animals (LAB)

6 hearts were harvested from surplus laboratory specific pathogen free Topig Norsvin pigs, mean weight 86 kg ± 21.5. As the body weight of laboratory animals is normally lower than the weight of slaughterhouse animals, a weight difference between those groups was unavoidable. All personnel was trained according to the EU Directive 2010/63/EU on the protection of animals used for scientific purposes, with adequate knowledge on hygiene measures to prevent zoonosis.

Anesthesia. Animals received premedication (midazolam intramuscular 0.4 mg/kg, ketamine intramuscular 10 mg/kg and atropine intramuscular 0.05 mg/kg) followed by thiopental (intravenous 4 mg/kg) for induction of anesthesia after which the animals were intubated. Anesthesia was maintained by midazolam (intravenous 1 mg/kg/h) and nimblex (intravenous 0.6 mg/kg/h).

Heart harvesting. The thorax and pericardium were opened. The pulmonary artery was dissected from the aorta and the root needle was inserted. 25,000 IU/L heparin was administered and after clamping the aorta, 1 L modified St Thomas cardioplegia crystalloid solution (added: Heparin 5000 IU/L) was administered with a target pressure of 100 mmHg and topical cooling was applied. The inferior caval vein was transected to unload the right atrium and the heart was excised.

Blood collection. As draining of 3.0–3.5 L of blood in the laboratory animals was impossible, the required amount of blood was collected from a SLH pig as mentioned before to assess the function of the heart on the PhysioHeart™ platform as described below.

Group difference. The difference in the animals per group were due to uncertainty of COMP settings in one heart. To secure statistical analysis another heart was added to the SLH group. After thorough analysis of the COMP perfusion data, we decided the heart could be used for analysis.

Cold oxygenated machine perfusion (COMP)

All hearts were stored with cold oxygenated machine perfusion (COMP) for 4 hours (Figure 1). The heart box system (XVIVO Perfusion, modified Kidney Assist transport box) consists of a reservoir, a pulsatile driven rotary pump (60BPM) and control unit. The heart box was assembled, primed with homemade modified Steen heart solution at 4°C (Supplementary Table 1), and deaired carefully. The temperature was regulated and




Figure 1. Schematic representation of the study design.

maintained around 8°C by placing ice around the reservoir. The machine was set to a pressure of 5 mmHg to maintain low flow before the heart was connected in the slaughterhouse.

The aorta was cannulated and de-aired carefully while the heart was submerged in preservation solution. Perfusion pressure was set to 20–25 mmHg with a target flow of 100–200 mL/min; all hearts met these criteria and were transported to our laboratory. The hearts were perfused for 4 hours including preparations for the PhysioHeart™ model.

Working heart model - PhysioHeart™

The PhysioHeart™ platform enables loading of the left side of the heart, which allows for functional analysis of left ventricular function. The platform can thus be seen as a surrogate to transplantation, which is more expensive and complex to perform in an experimental setting. Reperfusion of the heart on the PhysioHeart™ platform has been previously described.¹⁰ The system contains a pre- and afterload chamber and has the possibility to perfuse both in Langendorff and working mode. The system was primed with autologous heparinized diluted whole blood in the SLH group and with non-autologous heparinized diluted whole blood in the LAB group (diluted with modified Krebs Henseleit Buffer) to obtain 5 L final volume with $26\% \pm 2.41$ hematocrit (Supplementary Table 2). The blood was kept at 38°C and physiological blood gas parameters were maintained using a gas mixture of oxygen and carbon dioxide with an FiO₂ of 25%. The perfusate was supplemented with glucose and insulin to maintain physiological levels.

The left atrium, aorta and pulmonary artery were cannulated and connected to the system, while both caval veins and azygotic vein were closed with a suture. The hearts were then perfused in Langendorff mode with a target aortic pressure of 70–80 mmHg. Shortly after restoring blood flow to the myocardium, irregular contractions were observed, and the hearts were defibrillated with 30–50 J. If heartrate was less than 80 BPM, ventricular pacing using temporary pace leads with a target rate of 100 BMP (Medtronic Inc., Minneapolis, MN, USA) was initiated. A recovery perfusion

interval of 60 minutes in Langendorff mode was performed prior to switching to working mode by applying preload to the heart by filling of the left atrium and left ventricle thereby enabling cardiac output (CO). The left atrial pressure (LAP) was kept between 10 and 20 mmHg and the afterload was adjusted to create a physiological setting of CO and mean aortic pressure (MAP).⁴ Inclusion criteria for further analysis were a CO >3.0 L with a LAP <20 mmHg despite dobutamine administration. If these criteria were met, hearts were maintained in working mode for another 3 hours and function was assessed at specified time points (90, 120, 180, and 240 minutes).

Data collection

Functional data. Functional data was recorded for 10 seconds with 1000 Hz to create 10,000 data points per time point. Directly recorded data included: time in seconds, cardiac output (CO), coronary flow (CF), left ventricular pressure (LVP), left atrial pressure (LAP) and mean aortic pressure (MAP). With LVP and time, maximum developed pressure (dP/dT max) and minimum developed pressure (dP/dT min) were calculated. To correct for the difference in heart weights, CO and CF were divided by the heart weight, which resulted in cardiac index (CI) and coronary flow index (CFI). Survival was defined by a CO >3 L with a MAP >60 mmHg and a LAP <15 mmHg.

Bloodgas analysis. Bloodgas analysis was performed with a VetScan iSTAT 1 bloodgas analyzer (Abaxis, Union City, CA, USA) with CG8 + cartridges (Abbott, Chicago, IL, USA). Analysis was performed during hypothermic perfusion (5, 90, and 240 minutes), after priming of the PhysioHeart™ setup and during normothermic perfusion (10, 30, 60, 90, 120, 180, and 240 minutes)

Perfusate samples. Perfusate samples were collected during COMP and normothermic testing in Lithium-Heparin (BD Vacutainer plastic Lithium-Heparin tube) and Serum tubes (BD vacutainer plastic serum tube), centrifuged at 2300G (2 × 10 minutes), aliquoted and stored at –80°C. Samples were collected during hypothermic perfusion (5, 90, and 240 minutes) and

normothermic perfusion (baseline, 30, 60, 90, 120, 180, and 240 minutes). As weight could impact the results, variables analyzed were normalized to the mean weight of 450 g of all hearts included (variable was divided by the weight of the heart and multiplied by the mean weight of 450 g).

Table 1. Baseline and weight gain.

	SLH	LAB	p-value
Number	7	6	
Female (%)	86	67	0.052
Warm ischemic time (sec)			
Mean (SD)	192 (\pm 34)	NA	
Median (IQR)	198 (150–217)	NA	
Harvesting weight (g)			
Mean (SD)	526 (\pm 54)	371 (\pm 93)	0.003
Median (IQR)	514 (490–550)	399 (262–442)	0.004
Preservation weight gain (g)			
Mean (SD)	54 (\pm 43)	17 (\pm 52)	0.193
Median (IQR)	71 (6–79)	15 (-27–60)	0.153
Testing weight gain (g)			
Mean (SD)	52 (\pm 34)	98 (\pm 46)	0.065
Median (IQR)	54 (16–85)	107 (59–125)	0.063

Statistical analysis

Data analysis was performed using SPSS (IBM SPSS Statistics 26.0.0.1 (32bit)). Data is presented in mean \pm standard deviation (SD) and median with interquartile range, as applicable. Differences between groups were calculated either with student t-test or Mann-Whitney U depending on the data distribution. A mixed ANOVA was used to analyze the impact of time and a Kaplan-Meier curve with a log rank test was performed to detect differences in survival. A p-value below 0.05 was found significant. Correlations were determined by a Pearson correlation coefficient.

Results

Warm ischemic time, baseline heart weight and Langendorff data

The warm ischemic time for the LAB group was non-existent, as hearts were harvested in a controlled beating state, as it is the case in a DBD setting. The warm ischemic time of the SHL group it was 192 (\pm 34) seconds. Harvesting weight was lower in the LAB group (LAB 371 g (\pm 93) vs the SLH group 526 g (\pm 54), $p = 0.003$). Weight gain was similar in both groups during COMP and normothermic testing (Table 1).




Figure 2. Perfusate samples during cold oxygenated machine perfusion. SLH: slaughterhouse group; LAB: laboratory animal group; * ≤ 0.05 ; ** ≤ 0.001 .




Figure 3. Correlation plot for lactate at 240 minutes of cold oxygenated machine perfusion and cardiac index at 240 minutes of functional testing.

Table 2. Langendorff data.

	SLH	LAB	p-value
T02			
MAP (mmHg)	75 (± 5)	77 (± 3)	0.390
CF (L/min)	0.8 (± 0.1)	0.6 (± 0.1)	0.022
CFI (ml/min/g)	1.54 (± 0.32)	1.64 (± 0.43)	0.671
T30			
MAP (mmHg)	75 (± 5)	77 (± 5)	0.477
CF (L/min)	0.8 (± 0.2)	0.7 (± 0.1)	0.030
CFI (ml/min/g)	1.62 (± 0.35)	1.74 (± 0.49)	0.591
T59			
MAP (mmHg)	66 (± 9)	69 (± 7)	0.650
CF (L/min)	1.1 (± 0.1)	0.9 (± 0.1)	0.015
CFI (ml/min/g)	2.19 (± 0.33)	2.50 (± 0.57)	0.242

Perfusate samples during COMP

Lactate remained significantly lower in LAB hearts at 240 minutes ($1.0 (\pm 0.2)$ vs $4.9 (\pm 0.8)$ $p = 0.000$). This was also observed for lactate dehydrogenase and ammonia, whereas troponin levels did not differ. Oxygen levels decreased over time, with a steeper decline in the SLH hearts (Figure 2).

Lactate levels at 240 minutes of COMP correlated with the CI at 240 minutes of functional testing ($r = -0.868$, $p = < 0.000$) (Figure 3).

Baseline blood gas PhysioHeart™

No statistical differences were observed in most blood gas parameters at the start of normothermic perfusion on the PhysioHeart™ platform besides pO_2 (LAB 139 (± 31) vs SLH 176 (± 25), $p = 0.040$) and potassium (LAB

$4.4 (\pm 0.3)$ vs SLH $4.0 (\pm 0.3)$, $p = 0.030$) (Supplementary Table 2). One could argue whether these difference at the start of reperfusion are clinically relevant.

Langendorff mode. No significant difference was seen in perfusion pressure during Langendorff perfusion. The difference found in CF seems weight dependent as the CFI is not significantly different (Table 2).

Perfusate samples during physioHeart™

No differences were seen in lactate levels during normothermic testing, although lactate dehydrogenase was higher in the SLH group. Ammonia showed a trend towards higher levels in the LAB group, whereas no difference was found in troponin levels between groups (Figure 4).

Functional survival

Functional survival was similar in both groups, but LAB hearts showed significant better CI, CFI, LAP and developed pressures (Figure 5). Mixed ANOVA showed a significant time effect in functional decline and significant differences between groups for CO, CI, developed pressures, LAP, CFI and MAP (Supplementary Table 3).

Discussion

In this study, we compared 4 hours of COMP of SLH and LAB hearts to validate the use of SLH hearts for *ex-situ* machine perfusion studies. Functional survival in both groups was 100% during normothermic testing, but although LAB hearts had a superior function and required less dobutamine, both groups displayed a similar decline in functional parameters over time. The SLH hearts can therefore be used as an alternative to LAB hearts for *ex-situ* heart perfusion studies, as they respond in a similar way to the exposure of the perfusion setup with comparable survival and trends of function over time. The difference in function and differential requirement of dobutamine should be considered when analyzing data from SLH heart experiments.

Earlier studies have shown successful preservation and survival of porcine hearts from experimental animals for up to 24 hours using normothermic perfusion¹¹ and COMP.¹²⁻¹⁴ However, this paper is the first that describes the survival of SLH hearts for 4 hours. SLH hearts accessibility is unlimited, and they concomitantly reduce the cost of experiments. SLH hearts are harvested




Figure 4. Perfusate samples during normothermic perfusion and testing. SLH: slaughterhouse group; LAB: laboratory animal group; * ≤ 0.05 .

without buying and housing laboratory animals, the use of operating room, medication and extra personnel needed during the harvesting procedure.

Preserving hearts using COMP instead of SCS, which is still the golden standard in heart transplantation, is based on our hypothesis that COMP is a better preservation method. Key factors that support this hypothesis are a reduction in metabolism due to hypothermia, the washout or dilution of waste products and supplying oxygen, which has also been demonstrated in other hypothermically perfused organs.¹⁵⁻¹⁷ The washout of lactate, lactate dehydrogenase and ammonia were higher in the SLH hearts, whereas rise in troponin level was similar. These findings could be a result of the harvesting procedure in the SLH hearts with a period of warm ischemia, resulting in suboptimal preservation of energy stores prior to COMP. Although a higher level of damage markers was seen during COMP, all hearts showed clearance of lactate during normothermic perfusion and acceptable functional survival for 4 hours.

Interestingly, lactate levels at 240 minutes of COMP correlated with CI after 240 minutes of functional testing (Figure 3). At the start of clinical normothermic *ex-situ* heart perfusion, lactate was seen as a possible predictor of outcome after transplantation¹⁸ and was used as a marker to accept or deny hearts during normothermic perfusion.¹⁹ More recent research has shown that lactate during normothermic perfusion is

not very strongly correlated with functional outcomes.^{20,21} As a marker of adequate perfusion, the decrease in lactate provides sufficient information, for that reason no venous lactate was measured. Lactate as predictor for cardiac function during COMP is not well documented and needs further research.

Although SLH hearts showed functional survival for 4 hours, their function was inferior to LAB hearts. The superior cardiac function of LAB hearts can be attributed to the fact that these hearts were harvested in an optimal setting, resulting in the absence of warm ischemia. The length of the warm ischemic interval followed by reperfusion is inversely associated with cardiac function.²² Other important factors contributing to the functional difference between LAB and SLH animals are the electrical stunning of SLH animals that potentially contributes to more myocardial injury, the more extensive stress of these animals prior to harvesting and the finding that the SLH hearts are more hypertrophic compared to hearts from LAB animals. Hypertrophy was not precisely measured in this setup but was based on the inspection of axial coupes at the end of the experiment. The main factor that influences the hypertrophy is the fact that SLH animals are bred for human consumption. Farmers aim to reach a weight of around 120 kg by the age of 6 months with specific feeding strategies to optimize growth, unfortunately the weight of the SLH animals could not be measured due to logistic reasons. This is




Figure 5. Functional outcome. SLH: slaughterhouse group; LAB: laboratory animal group; dP/dTmax: maximum developed pressure; dP/dTmin: minimum developed pressure; * ≤ 0.05 ; ** ≤ 0.001 .

in contrast with LAB animals, which are not fed to gain weight in a short timeframe. The difference in harvesting weight was suboptimal but unavoidable. Laboratory animals are normally used at a younger age and are not bred or fed to gain weight in a short timeframe, which make animals lighter at similar age. Besides, the housing and handling of 120 kg pigs is

uncommon in the laboratory setting. For this reason, we decided to use hearts from animals that are commonly used in our animal facility and normalized perfusate samples to the mean weight of the group as described in our methods.

Despite the inferior function of SLH hearts, the trend of cardiac function over time showed to be

comparable for both groups. This reflects a similar response to reperfusion and normothermic testing. On the one hand this is a positive outcome, as this highlights again that SLH hearts can be used for *ex-situ* perfusion studies to optimize machine perfusion or test therapeutics. On the other hand, it highlights that potentially the reperfusion and normothermic testing for all hearts can be optimized. Literature describes this decline in function over time,²³ with metabolic pathways as potential target.

Another limiting factor in this study is the use of non-autologous blood during normothermic testing in the LAB hearts. Although we attempted to harvest the required amount of blood, we were not able to collect enough blood from the LAB animals. Since a minimal amount of full blood is needed for the normothermic testing phase, we have chosen to use non-autologous blood. The use of non-autologous blood in prior non-transfused pigs normally does not lead to hemolytic transfusion, but can have an impact on endothelial function, coagulation and inflammation.²⁴ Together this can impact cardiac function, although the extent of which is unknown. It could be that the decrease in cardiac function over time is partly due to the use of non-autologous blood, although similar declining trends were observed in the SLH hearts that were run with autologous blood. Although the use of non-autologous blood could have had an impact in the LAB group, cardiac function in the LAB group remained superior to the SLH hearts.

In conclusion, this study demonstrates that the SLH heart can be used as an alternative to LAB hearts for *ex-situ* perfusion studies. Although LAB hearts showed superior function at a lower dobutamine concentration, both LAB and SLH hearts showed similar functional survival and had a similar trend in functional decline during normothermic testing. The use of SLH hearts can reduce the need for laboratory animals and enables high volume access to porcine hearts for *ex-situ* perfusion studies in a cost-effective manner.


Declaration of conflictintg interests

The author(s) declared no potential conflicts of interest with respect to the research, authorship, and/or publication of this article.

Funding

The author(s) disclosed receipt of the following financial support for the research, authorship, and/or publication of this article: This work was supported by the Regenerative Medicine crossing borders (RegMed XB). This work was supported by European Research Council (ERC) under the EVICARE grant (number 725,229) to JS.

ORCID iD

Selma E Kaffka genaamd Dengler  <https://orcid.org/0000-0002-7407-8703>

Supplemental Material

Supplemental Material for this article is available online.

References

1. Messer S, Cernic S, Page A, et al. A 5-year single-center early experience of heart transplantation from donation after circulatory-determined death donors. *J Heart Lung Transplant* 2020; 39: 1463–1475. 2020/11/30. DOI: [10.1016/j.healun.2020.10.001](https://doi.org/10.1016/j.healun.2020.10.001)
2. Dhital K, Ludhani P, Scheuer S, et al. DCD donations and outcomes of heart transplantation: the Australian experience. *Indian J Thorac Cardiovasc Surg* 2020; 36: 224–232. 2020/10/17. DOI: [10.1007/s12055-020-00998-x](https://doi.org/10.1007/s12055-020-00998-x)
3. Prescott MJ, Lidster K. Improving quality of science through better animal welfare: the NC3Rs strategy. *Lab Anim (NY)* 2017; 46: 152–156. 2017/03/23. DOI: [10.1038/labani.1217](https://doi.org/10.1038/labani.1217)
4. Kappler B, Ledezma CA, van Tuijl S, et al. Investigating the physiology of normothermic ex vivo heart perfusion in an isolated slaughterhouse porcine model used for device testing and training. *BMC Cardiovasc Disord* 2019; 19: 254–2019. 2019/11/13. DOI: [10.1186/s12872-019-1242-9](https://doi.org/10.1186/s12872-019-1242-9)
5. Goumon S, Faucitano L. Influence of loading handling and facilities on the subsequent response to pre-slaughter stress in pigs. *Livestock Science* 2017; 200: 6–13. DOI: [10.1016/j.livsci.2017.03.021](https://doi.org/10.1016/j.livsci.2017.03.021)
6. Martínez-Miró S, Tecles F, Ramón M, et al. Causes, consequences and biomarkers of stress in swine: an update. *BMC Vet Res* 2016; 12: 171. 171 2016/08/21. DOI: [10.1186/s12917-016-0791-8](https://doi.org/10.1186/s12917-016-0791-8)
7. EU. Council regulation (EC) 1099/2009 on the protection of animals at the time of killing 2009: 1–33.
8. Gaasch WH, Zile MR, Hoshino PK, et al. Tolerance of the hypertrophic heart to ischemia. studies in compensated and failing dog hearts with pressure overload hypertrophy. *Circulation* 1990; 81: 1644–1653. 1990/05/01. DOI: [10.1161/01.cir.81.5.1644](https://doi.org/10.1161/01.cir.81.5.1644)
9. Copeland H, Hayanga JWA, Neyrinck A, et al. Donor heart and lung procurement: a consensus statement. *J Heart Lung Transplant* 2020; 39: 501–517. 2020/06/07. DOI: [10.1016/j.healun.2020.03.020](https://doi.org/10.1016/j.healun.2020.03.020)
10. de Hart J, de Weger A, van Tuijl S, et al. An ex vivo platform to simulate cardiac physiology: a new dimension for therapy development and assessment. *Int J Artif Organs* 2011; 34: 495–505. 2011/07/05. DOI: [10.5301/ijao.2011.8456](https://doi.org/10.5301/ijao.2011.8456)
11. Tchouta L, Drake D, Hoenerhoff M, et al. Twenty-four-hour normothermic perfusion of isolated ex vivo hearts using plasma exchange. *J Thorac Cardiovasc Surg* 2022; 164: 128–138. 2021/01/25. DOI: [10.1016/j.jtcvs.2020.11.158](https://doi.org/10.1016/j.jtcvs.2020.11.158)
12. Qin G, Sjöberg T, Liao Q, et al. Intact endothelial and contractile function of coronary artery after 8 hours of

- heart preservation. *Scand Cardiovasc J* 2016; 50: 362–366. 2016/07/16. DOI: [10.1080/14017431.2016.1213876](https://doi.org/10.1080/14017431.2016.1213876)
13. Steen S, Paskevicius A, Liao Q, et al. Safe orthotopic transplantation of hearts harvested 24 hours after brain death and preserved for 24 hours. *Scand Cardiovasc J* 2016; 50: 193–200. 2016/02/18. DOI: [10.3109/14017431.2016.1154598](https://doi.org/10.3109/14017431.2016.1154598)
 14. Michel SG, La Muraglia GM 2nd, Madariaga ML, et al. Twelve-Hour Hypothermic Machine Perfusion for Donor Heart Preservation Leads to Improved Ultrastructural Characteristics Compared to Conventional Cold Storage. *Ann Transplant* 2015; 20: 461–468. 2015/08/12. DOI: [10.12659/aot.893784](https://doi.org/10.12659/aot.893784)
 15. Rivard AL, Gallegos R, Ogden IM, et al. Perfusion preservation of the donor heart: basic science to pre-clinical. *J Extra Corpor Technol* 2009; 41: 140–148. DOI: 2009/10/08
 16. Dutkowski P, Schlegel A, de Oliveira M, et al. HOPE for human liver grafts obtained from donors after cardiac death. *J Hepatol* 2014; 60: 765–772. 2013/12/04. DOI: [10.1016/j.jhep.2013.11.023](https://doi.org/10.1016/j.jhep.2013.11.023)
 17. Jochmans I, Brat A, Davies L, et al. Oxygenated versus standard cold perfusion preservation in kidney transplantation (COMPARE): a randomised, double-blind, paired, phase 3 trial. *Lancet (London, England)* 2020; 396: 1653–1662. 2020/11/23. DOI: [10.1016/s0140-6736\(20\)32411-9](https://doi.org/10.1016/s0140-6736(20)32411-9)
 18. Hamed A, Tsui S, Huber J, et al. A Serum lactate is a highly specific predictor of post cardiac transplant outcomes using the organ care system. *J Heart Lung Transplant* 2009; 28: S71.
 19. Ardehali A, Esmailian F, Deng M, et al. Ex-vivo perfusion of donor hearts for human heart transplantation (PROCEED II): a prospective, open-label, multicentre, randomised non-inferiority trial. *Lancet* 2015; 385: 2577–2584. 2015/04/19. DOI: [10.1016/s0140-6736\(15\)60261-6](https://doi.org/10.1016/s0140-6736(15)60261-6)
 20. White CW, Ambrose E, Müller A, et al. Assessment of donor heart viability during ex vivo heart perfusion. *Can J Physiol Pharmacol* 2015; 93: 893–901. 2015/09/01. DOI: [10.1139/cjpp-2014-0474](https://doi.org/10.1139/cjpp-2014-0474)
 21. Messer SJ, Axell RG, Colah S, et al. Functional assessment and transplantation of the donor heart after circulatory death. *J Heart Lung Transplant* 2016; 35: 1443–1452. 2016/12/06. DOI: [10.1016/j.healun.2016.07.004](https://doi.org/10.1016/j.healun.2016.07.004)
 22. Frank A, Bonney M, Bonney S, et al. Myocardial ischemia reperfusion injury: from basic science to clinical bedside. *Semin Cardiothorac Vasc Anesth* 2012; 16: 123–132. 2012/03/01. DOI: [10.1177/1089253211436350](https://doi.org/10.1177/1089253211436350)
 23. Hatami S, White CW, Shan S, et al. Myocardial functional decline during prolonged ex situ heart perfusion. *Ann Thorac Surg* 2019; 108: 499–507. 2019/03/16. DOI: [10.1016/j.athoracsur.2019.01.076](https://doi.org/10.1016/j.athoracsur.2019.01.076)
 24. Smith DM, Newhouse M, Naziruddin B, et al. Blood groups and transfusions in pigs. *Xenotransplantation* 2006; 13: 186–194. 2006/06/08. DOI: [10.1111/j.1399-3089.2006.00299.x](https://doi.org/10.1111/j.1399-3089.2006.00299.x)

Appendix

Nomenclature

BPM	Beats per minute
CF	Coronary flow
CFI	Coronary flow index
CI	Cardiac index
CO	Cardiac output
COMP	Cold oxygenated machine perfusion
dP/dT max	Maximum developed pressure
dP/dT min	Minimum developed pressure
IRI	Ischemia reperfusion injury
LAB	Laboratory
LAP	Left atrial pressure
LVP	Left ventricular pressure
MAP	Mean aortic pressure
SLH	Slaughterhouse.